

Application	DNA Sequencing										Metagenomics		Targeted Sequencing	
	De Novo Assembly - HFI Reads	De Novo Assembly - for Low DNA Input	De Novo Assembly and Variant Detection - for Ultra-Low DNA Input	Microbial De Novo Assembly	Variant Detection	Structural Variation Detection	Iso-Seq Method	Single-Cell Iso-Seq	Full-length 16S rRNA Sequencing	Shotgun Metagenomic Profiling or Assembly	Amplicon Sequencing	No-Amp Targeted Sequencing		
<b>Experimental Design</b>														
<b>For Sequel II Systems Applications-specific Insert Size and Read Numbers, see the Purchasing Guide</b>														
<b>With 1 SMRT Cell you can:</b>	Produce reference quality assemblies for genomes up to 2 Gb	Produce reference quality assemblies for genomes up to 1 Gb. Multiplex up to 2 small genomes on the Sequel II System	Produce reference quality assemblies for genomes up to 500 Mb	Sequence up to 48 microbes	With 2 SMRT Cells, SM, Cell SVs, iCNVs, and SVs in a 3 Gb genome	Call SVs for up to 2 samples with ~9 Gb genomes	Characterize alternative splicing/transcribe a genome with full length transcripts	Characterize alternative splicing with full length transcripts up to 3M full length reads with cell barcode and UMI information	Multiplex up to 96 samples to provide strain level resolution	Generate near-complete assemblies of high-complexity samples (e.g. gut microbiome)	Sequence 384 barcoded amplicons	Sequence 5 targeted regions in a multiplex of 20 samples		
<b>Minimum Recommended Coverage</b>	10- to 15-fold HFI read coverage per haplotype	10- to 15-fold HFI read coverage per haplotype	>30-fold HFI read coverage per diploid sample for de novo assembly >15-fold HFI read coverage for human variant detection	≥30-fold Unique Molecular Coverage (UMC) per microbial genome	≥15-fold HFI read coverage of a human genome	5- to 20-fold UMC per diploid sample depending on study goals	One human transcriptome per SMRT Cell BM	1,000 unique reads/ single cell for 3000 cells 10,000 unique reads/ single cell for 300 cells	8,000 HFI reads / 16S sample	See Best practices guide	≥150-fold HFI read coverage per target locus for variant detection 6,000-fold HFI read coverage per target locus for minor variant detection (1% sensitivity)	≥100-fold HFI read coverage per target locus		
<b>Library Insert Size</b>	15 - 20 kb	~15 kb	10 - 12 kb	10 - 15 kb	15 - 20 kb	≥15 kb	<2 kb to >3 kb	<2 kb to >3 kb	1 - 2 kb	10 kb	500 bp - 15 kb	4-8 kb or larger		
<b>Sample Preparation</b>														
<b>Procedure and Checklist Reference</b>	<a href="#">Preparing HFI SMRTbell Libraries Using PacBio Sequencing Template Prep Kit v2.0</a>	<a href="#">Preparing HFI Libraries from Low DNA Input Using SMRTbell Express Template Prep Kit v2.0</a>	<a href="#">Preparing HFI SMRTbell Libraries from Ultra-Low DNA Input</a>	<a href="#">Preparing Microbial SMRTbell Libraries Using PacBio Sequencing Template Prep Kit v2.0</a>	<a href="#">Preparing HFI SMRTbell Libraries Using PacBio Sequencing Template Prep Kit v2.0 for Structural Variant Detection</a>	<a href="#">Preparing gRNA Libraries Using the SMRTbell Express Template Preparation Kit</a>	<a href="#">Iso-Seq Express Template Preparation for Sequel and Sequel II Systems</a>	<a href="#">Preparing Single-Cell Iso-Seq Libraries Using SMRTbell Express Template Prep Kit v2.0</a>	<a href="#">Amplification of Full-Length 16S Gene with Barcoded Primers for Multiplexed SMRTbell Library Preparation and Sequencing</a>	<a href="#">Preparing 16S Libraries Using SMRTbell Express Template Prep Kit v2.0 for Metagenomic Strain Sequencing</a>	<a href="#">Preparing SMRTbell Libraries Using PacBio Sequencing Template Prep Kit v2.0 for Targeted Sequencing</a>	<a href="#">Preparing SMRTbell Libraries Using the PCR-Free method</a>		
<b>Minimum Input Amount</b>	15 µg	>400 ng per 1 Gb genome size (single-sample) >800 ng per 600 Mb genome size (2-pile)	5 - 20 ng per 500 Mb genome size	1 µg per microbe	15 µg	3 µg	300 ng total RNA for 1st Strand cDNA Synthesis	>100 ng cDNA AFTER reamplification	500 ng - 1 µg	1.5 µg	250-500 ng for 250-1000 bp 500-1000 ng for 1-3 kb 1000-2000 ng for 3-10 kb 1500-3000 ng for 15kb	5 to 10 µg (represented by either a single sample or the total of multiple samples that will be multiplexed)		
<b>Recommended Pacific Biosciences Template Prep Kit</b>	Express TPK 2.0 + SMRTbell Enzyme Cleanup kit + Sequencing primer v2	Express TPK 2.0 + Sequencing primer v4	Express TPK 2.0 + Sequencing primer v4	Express TPK 2.0 + Sequencing primer v4	Express TPK 2.0 + SMRTbell Enzyme Cleanup kit + Sequencing primer v2	Express TPK 2.0 + Sequencing primer v2	Express TPK 2.0 + Sequencing primer v4	Express TPK 2.0 + Sequencing primer v4	Express TPK 2.0 + Sequencing primer v4	Express TPK 2.0 + Sequencing primer v2	Express TPK 2.0 + Sequencing primer v4	No-Amp Accessory Kit		
<b>Multiplexing/SMRT Cell</b>	N/A	N/A	N/A	Up to 48 microbes / SMRT Cell 1M	N/A	Up to 2 human samples/ SMRT Cell 1M	The protocol supports up to 12 barcodes available.	Detects cell barcodes and UMIs	Up to 96 samples/ SMRT Cell BM	Profile one community/ SMRT Cell BM	Up to 1,000+ samples/ SMRT Cell BM or SMRT Cell 1M	Up to 48 samples/ SMRT Cell BM		
<b>SMRT Sequencing with the Sequel II and Sequel IIe Systems: Loading and Pre-Extension Recommendations</b>														
<b>Sequencing Performance and Read Characteristics</b>	<b>Pacific Biosciences Sequel II Binding Kit</b>	2.0	2.0	2.0	2.0	2.0	2.0	2.0 / 2.1*	2.0 / 2.1*	2.1	2.0	2.1 - 500bp - 3,000 bp 2.0 - 20,000 bp	2.0	
	<b>Sequel II Sequencing Plate</b>	2.0	2.0	2.0	2.0	2.0	2.0	2.0	2.0	2.0	2.0	2.0	2.0	
	<b>Sequencing Mode</b>	CCS	CCS	CCS	CLR	CCS	CLR	CCS	CCS	CCS	CCS	CCS	CCS	
	<b>Movie Collection Time</b>	30 h	15 / 30 h	30 h	15 h	30 h	15 h	24 h	24 h	10 h	30 h	Insert Size Dependent	≥10 h (30 h for repeat expansion targets)	
	<b>Notes</b>	*Use Binding Kit 2.1 as default and 2.0 for long transcripts												
<b>Sequencing Performance and Read Characteristics</b>	<b>SMRT Sequencing with the Sequel IIe System: Loading and Pre-Extension Recommendations</b>													
	<b>Pacific Biosciences Sequel IIe Binding Kit</b>	3.0	3.0	NA	3.0	3.0	3.0	3.0	3.0	3.0	3.0	3.0	3.0	
	<b>Sequel IIe Sequencing Plate</b>	3.0	3.0	NA	3.0	3.0	3.0	3.0	3.0	3.0	3.0	3.0	3.0	
	<b>Sequencing Mode</b>	CCS	CCS	NA	CLR	CCS	CLR	CCS	CCS	CCS	CCS	CCS	CCS	
	<b>Movie Collection Time</b>	20 h	20 h	NA	10 h	20 h	10 h	20 h	20 h	10 h	20 h	Insert Size Dependent	≥10 h (20 h for repeat expansion targets)	
<b>Notes</b>	NA													
<b>Data Analysis Tools Available Through SMRT Link, PacBio DevNet and Other Compatible Software Tools</b>														
<b>Other Compatible Software Tools for Sequencing Data</b>	<b>SMRT Analysis GUI Applications</b>	Genome Assembly powered by IPA	Genome Assembly powered by IPA	Mark PCR Duplicates, Trim gRNA Amplification Adapter, Structural Variant Calling, Genome Assembly powered by IPA, imz_gbamarkdup.pbay	Demultiplex Barcodes followed by Microbial Assembly analysis	CCS with Mapping	Structural Variant Calling	Iso-Seq	Iso-Seq	CCS	CCS	CCS or Long Amplicon Analysis	Demultiplex Barcodes followed by CCS with Mapping analysis	
	<b>PacBio DevNet Tools</b>													
	<b>Other Analysis Tools</b>	hifiasm, HiCanu	hifiasm, HiCanu	hifiasm, HiCanu	Recommended: Google DeepVariant								<a href="#">GitHub Tools</a>	<a href="#">PacBio Amplicon Analysis (aba)</a>

Read lengths, reads/data per SMRT Cell and other sequencing performance results vary based on sample quality/type and insert size

[PacBio Glossary of Terms](#)

For Research Use Only. Not for use in diagnostic procedures. © Copyright 2020, Pacific Biosciences of California, Inc. All rights reserved. Information in this document is subject to change without notice. Pacific Biosciences assumes no responsibility for any errors or omissions in this document. Certain notices, terms, conditions and/or use restrictions may pertain to your use of Pacific Biosciences products and/or third party products. Please refer to the applicable Pacific Biosciences Terms and Conditions of Sale and to the applicable license terms at <https://www.pacb.com/legal-and-trademarks/terms-and-conditions-of-sale/>. Pacific Biosciences, the Pacific Biosciences logo, PacBio, SMRT, SMRTbell, Iso-Seq and Sequel are trademarks of Pacific Biosciences. BluePippin and SageELF are trademarks of Sage Science, Inc. NGS-seq and NGSengine are trademarks of Genix. Femo Pulse and Fragment Analyzer are trademarks of Agilent Technologies, Inc. All other trademarks are the sole property of their respective owners. Pacific Biosciences does not sell a kit for carrying out the overall No-Amp Targeted Sequencing method. Use of this method may require rights to third-party owned intellectual property. PN 101451-300 version 07 (Dec 2020)

