



Aria-MBI 85um nozzle 103008

Cytometer:	FACSAriaII	User:	Craig
Cytometer Name:	FACSAriaII	Institution:	ICBR Staff
Serial Number:	P46900101	Software:	BD FACSDiva 8.0.1
Input Device:	Manual	Date:	11/30/2018 10:58:39 AM
Sheath Pressure:	30.00		
Nozzle Size:	85		
Window Extension:	2.00		


Red Laser (640nm)




Green Laser
(532nm)




Blue Laser
(488nm)
FSC



UV Laser
(355nm)



Violet Laser
(405nm)



Comments:

Original BD configuration @ 45 PSI does not work so well. Changed to 30psi.
NAB, 10/30/2008